Managing High Ketones Pump Therapy

When Blood Glucose >14mmols/l it is essential to check ketone level.

Always have a Pen Device and Penfil Insulin as back up and part of your Emergency Pack.

If you are not sure what to do please give us a call early to avoid an admission to hospital.

Daytime (8am-6pm) 0151 252 5766 Out of Hours 0151 228 4811 Ask for the Diabetes Nurse on call service

Managing High Ketones Pump Therapy

When Blood Glucose >14mmols/l it is essential to check ketone level

Blood Ketone Guide		
Less than 0.6mmol/L	0.6-1.5mmol/L	More than 1.5mmol/L
Give a normal correction dose	Increase correction dose by 50%	Double correction dose

Novorapid, Humalog or Apidra is to be used for correction doses.

The correction dose must not exceed more than 20% of your Total Daily Insulin Dose.

REMEMBER!

- If blood glucose 14mmol or above with ketones 0.6mmol/l or above additional insulin MUST be given by a Pen Device NOT your Pump.
- CHANGE the insulin, cannula & infusion set
- START a 30%-60% temporary basal increase
- RECHECK in 1 hour

If vomiting or change in breathing attend A&E



NHS Foundation Trust